Material Safety Data Sheet

Date prepared: May 28, 2012 Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name: Sulfathiazole
CAS Reg. Number: 72-14-0
Catalog Number: 1ST4007
Manufacturer's Name: A Chemtek, Inc.

100 Barber Avenue

Worcester, MA 01606, USA

Telephone Number: 508-471-4121 **Fax Number:** 508-845-9201

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

Emergency Overview

OSHA Hazards Target Organ Effect, Irritant

Target OrgansBlood, Kidney **GHS Classification**

Acute toxicity, Oral (Category 5) Skin irritation (Category 2) Eye irritation (Category 2A)

Specific target organ toxicity - single exposure (Category 3)

GHS Label elements, including precautionary statements



Pictogram
Signal word
Warning

Hazard statement(s)

H303 May be harmful if swallowed.
 H315 Causes skin irritation.
 H319 Causes serious eye irritation.
 H335 May cause respiratory irritation.

Precautionary statement(s)

P261 Avoid breathing dust.

P264 Wash skin thoroughly after handling.

P271 Use only outdoors or in a well-ventilated area.

P280 Wear protective gloves/ eye protection/ face protection.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact

lenses, if present and easy to do. Continue rinsing.

P312 Call a POISON CENTER or doctor/ physician if you feel unwell.
P321 Specific treatment (see supplemental first aid instructions on this label).

P332 + P313 If skin irritation occurs: Get medical advice/ attention.
P337 + P313 If eye irritation persists: Get medical advice/ attention.
P362 Take off contaminated clothing and wash before reuse.
P403 + P233 Store in a well-ventilated place. Keep container tightly closed.

P405 Store locked up.

P501 Dispose of contents/ container to an approved waste disposal plant.

HMIS Classification

Health hazard: 2

Chronic Health Hazard: *

Flammability: 0 Physical hazards: 0

NFPA Rating

Health hazard: 2

Fire: 0

Reactivity Hazard: 0

Material Safety Data Sheet

Date prepared: May 28, 2012 Date revised: N/A (first edition)

Section 3 - Composition/Information on Ingredients

Formula: $C_9H_9N_3O_2S_2$ Molecular Weight: 255.32g/molCAS-No. 72-14-0

Section 4 - First Aid Measures

If inhaled

If breathed in, move person into fresh air. If not breathing give artificial respiration.

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Flush eyes with water as a precaution.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

Section 5 - Fire-Fighting Measures

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self contained breathing apparatus for fire fighting if necessary.

Hazardous combustion products

Hazardous decomposition products formed under fire conditions. - Carbon oxides, nitrogen oxides (NOx), Sulphur oxides

Section 6 - Accidental Release Measures

Personal precautions Avoid dust formation.

Environmental precautions Do not let product enter drains.

Methods for cleaning up Sweep up and shovel. Keep in suitable, closed containers for disposal.

Section 7 - Handling and Storage

Handling Avoid contact with skin and eyes. Avoid formation of dust and aerosols.

Provide appropriate exhaust ventilation at places where dust is formed. Normal measures for

preventive fire protection.

Storage Keep container tightly closed in a dry and well-ventilated place.

Section 8 - Exposure Controls/Personal Protection

Contains no substances with occupational exposure limit values.

Personal protective equipment

Eye/face protection

Safety glasses with side-shields conforming to EN166 Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin protection

Impervious clothing, The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Body Protection

The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Section 9 - Physical and Chemical Properties

Appearance:

Form Powder Colour White

Safety data

Material Safety Data Sheet

Date prepared: May 28, 2012 Date revised: N/A (first edition)

pH: not available

Melting point: Melting point/range: 200 - 202 °C (392 - 396 °F) - lit.

Boiling point: not available
Flash point: not available
Ignition temperature: not available
Lower explosion limit: not available
Upper explosion limit: not available
Water solubility: not available
Other safety information not available

Section 10 - Stability and Reactivity

Storage stability Stable under recommended storage conditions

Materials to avoid Strong oxidizing agents

Hazardous decomposition products Hazardous decomposition products formed under fire conditions. -

Carbon oxides, nitrogen oxides (NOx), Sulphur oxides Other decomposition products - no data available

Section 11 - Toxicological Information

Acute toxicity

LD50 Oral – mouse – 4.500 mg/kg

Irritation and corrosion:not availableSensitization:not availableRespiratory or skin sensitization:not available

Chronic exposure

No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC, ACGIH, NTP, or OSHA.

Potential health effects

Inhalation: May be harmful if inhaled. May cause respiratory tract irritation.

Skin: May be harmful if absorbed through skin. May cause skin irritation.

Eyes: May cause eye irritation.
Ingestion: May be harmful if swallowed.

Section 12 - Ecological Information

Toxicity not available

Elimination information (persistence and degradability)

not available

Ecotoxicity effects: not available **Further information on ecology:** not available

Section 13 - Disposal Considerations

Product Observe all federal, state, and local environmental regulations.

Contaminated packaging Dispose of as unused product.

Section 14- Transport Information

DOT (US):

IMDG:

Not dangerous goods

Not dangerous goods

Not dangerous goods

Not dangerous goods

Section 15 - Regulatory Information

OSHA Hazards

Target Organ Effect, Irritant

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the

Material Safety Data Sheet

Date prepared: May 28, 2012 Date revised: N/A (first edition)

threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

Acute Health Hazard, Chronic Health Hazard

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Sulfathiazole CAS-No.72-14-0

New Jersey Right To Know Components

Sulfathiazole CAS-No.72-14-0 California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.