Material Safety Data Sheet

Date prepared: Sep 5, 2015 Date revised: N/A (first edition)

Section 1 - General Information

Chemical Name:LinaloolCAS Reg. Number:78-70-6Catalog Number:1ST3303Manufacturer's Name:A Chemt

A Chemtek, Inc. 100 Barber Avenue

Worcester, MA 01606, USA

Telephone Number: 508-471-4121 **Fax Number:** 508-845-9201

For R & D use only

Section 2 - Hazardous Ingredients/Identity Information

Classification of the substance or mixture

Flammable liquids (Category 4), H227 Skin irritation (Category 2), H315 Eye irritation (Category 2A), H319 Acute aquatic toxicity (Category 3), H402



Pictogram
Signal word
Warning

Hazard statement(s)

H227 Combustible liquid.
H315 Causes skin irritation.
H319 Causes serious eye irritation.
H402 Harmful to aquatic life.

Precautionary statement(s)

P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking.

P264 Wash skin thoroughly after handling. P273 Avoid release to the environment.

P280 Wear protective gloves/ eye protection/ face protection.
P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove

contact lenses, if present and easy to do. Continue rinsing.

P332 + P313 If skin irritation occurs: Get medical advice/ attention.
P337 + P313 If eye irritation persists: Get medical advice/ attention.
P362 Take off contaminated clothing and wash before reuse.

P370 + P378 In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish.

P403 + P235 Store in a well-ventilated place. Keep cool.

P501 Dispose of contents/ container to an approved waste disposal plant.

Section 3 - Composition/Information on Ingredients

Formula: C10H18O Molecular Weight: 154.25 g/mol CAS-No. 78-70-6

Section 4 - First Aid Measures

If inhaled

If breathed in, move person into fresh air. If not breathing give artificial respiration.

In case of skin contact

Wash off with soap and plenty of water.

In case of eye contact

Material Safety Data Sheet

Date prepared: Sep 5, 2015 Date revised: N/A (first edition)

Flush eyes with water as a precaution.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water.

Section 5 - Fire-Fighting Measures

Suitable extinguishing media

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

Special protective equipment for fire-fighters

Wear self-contained breathing apparatus for firefighting if necessary.

Section 6 - Accidental Release Measures

Personal precautions Wear respiratory protection. Avoid breathing vapours, mist or gas. Ensure

adequate ventilation. Evacuate personnel to safe areas. For personal protection see section 8.

Environmental precautions Prevent further leakage or spillage if safe to do so. Do not let product enter

drains. Discharge into the environment must be avoided.

Methods for cleaning up Soak up with inert absorbent material and dispose of as hazardous waste. Keep

in suitable, closed containers for disposal.

Section 7 - Handling and Storage

Handling Avoid contact with skin and eyes. Avoid inhalation of vapour or mist.

Storage Keep container tightly closed in a dry and well-ventilated place.

Containers which are opened must be carefully resealed and kept upright to prevent leakage.

Recommended storage temperature: 2 - 8 °C.

Section 8 - Exposure Controls/Personal Protection

Contains no substances with occupational exposure limit values

Personal protective equipment

Respiratory protection

Where risk assessment shows air-purifying respirators are appropriate use a full-face respirator with multipurpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Hand protection For prolonged or repeated contact use protective gloves

Eye protection Safety glasses

Hygiene measures General industrial hygiene practice

Section 9 - Physical and Chemical Properties

Appearance:

Form liquid Color colourless

Safety data

pH: not available
Melting point: not available
Boiling point: not available
Flash point: not available
Lower explosion limit: not available
Upper explosion limit: not available
Water solubility: not available

Section 10 - Stability and Reactivity

Storage stability Stable under recommended storage conditions

Material Safety Data Sheet

Date prepared: Sep 5, 2015 Date revised: N/A (first edition)

Materials to avoid Strong oxidizing agents, Strong bases

Hazardous decomposition productsHazardous decomposition products formed under fire conditions.

- Carbon oxides.

Section 11 - Toxicological Information

Acute toxicity: LD50 Oral - Rat - male and female - 2,790 mg/kg

LD50 Dermal - Rabbit - 5,610 mg/kg

Irritation and corrosion:not availableSensitization:not availableRespiratory or skin sensitization:not available

Chronic exposure

IARC: No component of this product present at levels greater than or equal to 0.1% is identified as

probable, possible or confirmed human carcinogen by IARC.

ACGIH: No component of this product present at levels greater than or equal to 0.1% is identified as a

carcinogen or potential carcinogen by ACGIH.

NTP: No component of this product present at levels greater than or equal to 0.1% is identified as a

known or anticipated carcinogen by NTP.

OSHA: No component of this product present at levels greater than or equal to 0.1% is identified as a

carcinogen or potential carcinogen by OSHA.

Potential health effects

Inhalation: no data available
Skin: no data available
Eyes: no data available
Ingestion: no data available

Section 12 - Ecological Information

Elimination information (persistence and degradability)

Toxicity to fish static test LC50 - Oncorhynchus mykiss (rainbow trout) - 27.8 mg/l - 96 h
Toxicity to daphnia and Immobilization EC50 - Daphnia magna (Water flea) - 59 mg/l - 48 h

other aquatic invertebrates

Toxicity to algae static test EC50 - Desmodesmus subspicatus (green algae) - 156.7 mg/l - 96 h

Ecotoxicity effects: not available **Further information on ecology:** not available

Section 13 - Disposal Considerations

Product Offer surplus and non-recyclable solutions to a licensed disposal company.

Contaminated packaging Dispose of as unused product.

Section 14- Transport Information

DOT (US):

NA-Number: 1993 Class: NONE Packing group: III Proper shipping name: Combustible liquid, n.o.s. (Linalool)

Reportable Quantity (RQ): Poison Inhalation Hazard: No

IMDG:

Not dangerous goods

IATA:

Not dangerous goods

Section 15 - Regulatory Information

Material Safety Data Sheet

Date prepared: Sep 5, 2015 Date revised: N/A (first edition)

SARA 302 Components

SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

SARA 313: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

Massachusetts Right To Know Components

No components are subject to the Massachusetts Right to Know Act.

Pennsylvania Right To Know Components

Linalool CAS-No. 78-70-6 Revision Date

New Jersey Right To Know Components

Linalool CAS-No. 78-70-6 Revision Date

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Section 16 - Other Information

The above information is believed to be correct but does not purport to be all-inclusive and shall be used only as a guide. A Chemtek shall not be held liable for any damage resulting from handling or from contact with the above product.

- 4 -